

Phospholysine phosphohistidine inorganic pyrophosphate phosphatase

Catalog No: #AP84925

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84925-1 50ug #AP84925-2 100ug #AP84925-3 1mg

Description

Product Name	Phospholysine phosphohistidine inorganic pyrophosphate phosphatase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Accession No.	Uniprot ID:Q5I0D5
Uniprot	Q5I0D5
GeneID	361663;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lhpp

Note: This product is for in vitro research use only