

Protein Lines homolog 1

Catalog No: #AP84926



Package Size: #AP84926-1 50ug #AP84926-2 100ug #AP84926-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein Lines homolog 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Wnt-signaling molecule Lines homolog 1,LINS,WINS1
Accession No.	Uniprot ID:Q8NG48
Uniprot	Q8NG48
GeneID	55180;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LINS1

Note: This product is for in vitro research use only