

Liver-expressed antimicrobial peptide 2

Catalog No: #AP84928



Package Size: #AP84928-1 50ug #AP84928-2 100ug #AP84928-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Liver-expressed antimicrobial peptide 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	LEAP-2
Accession No.	Uniprot ID:Q95JC3
Uniprot	Q95JC3
GeneID	282270;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LEAP2

Note: This product is for in vitro research use only