

Leptin receptor

Catalog No: #AP84929

Package Size: #AP84929-1 50ug #AP84929-2 100ug #AP84929-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Leptin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LEP-R,OB receptor,OB-R,CD295,Fa,Obr
Accession No.	Uniprot ID:Q62959
Uniprot	Q62959
GeneID	24536;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lepr

Note: This product is for in vitro research use only