

Induced myeloid leukemia cell differentiation protein Mcl-1 homolog

Catalog No: #AP84931

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84931-1 50ug #AP84931-2 100ug #AP84931-3 1mg

Description

Product Name	Induced myeloid leukemia cell differentiation protein Mcl-1 homolog
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Bcl-2-related protein EAT/mcl1
Accession No.	Uniprot ID:P97287
Uniprot	P97287
GeneID	17210;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mcl1

Note: This product is for in vitro research use only