

# Leukotriene C4 synthase

Catalog No: #AP84941



Package Size: #AP84941-1 50ug #AP84941-2 100ug #AP84941-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leukotriene C4 synthase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LTC4 synthase,Leukotriene-C(4) synthase
Accession No.	Uniprot ID:Q925U2
Uniprot	Q925U2
GeneID	114097;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Ltc4s

Note: This product is for in vitro research use only