

Ragulator complex protein LAMTOR5

Catalog No: #AP84942



Package Size: #AP84942-1 50ug #AP84942-2 100ug #AP84942-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ragulator complex protein LAMTOR5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	HBV X-interacting protein homolog,HBX-interacting protein homolog,Xip
Accession No.	Uniprot ID:Q9D1L9
Uniprot	Q9D1L9
GeneID	68576;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lamtor5

Note: This product is for in vitro research use only