

Lysyl oxidase homolog 4

Catalog No: #AP84948



Package Size: #AP84948-1 50ug #AP84948-2 100ug #AP84948-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lysyl oxidase homolog 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Lysyl oxidase-like protein 4,Lysyl oxidase-related protein C,LOXC
Accession No.	Uniprot ID:Q96JB6
Uniprot	Q96JB6
GeneID	84171;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LOXL4

Note: This product is for in vitro research use only