

# Krueppel-like factor 10

Catalog No: #AP84949

Package Size: #AP84949-1 50ug #AP84949-2 100ug #AP84949-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Krueppel-like factor 10
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	GDNF-inducible factor, Transcription factor GIF, mGIF, Transforming growth factor-beta-inducible early growth response protein 1, TGFB-inducible early growth response protein 1, Gdnfif, Tieg, Tieg1
Accession No.	Uniprot ID:O89091
Uniprot	O89091
GeneID	21847;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (PH8.0) gene\_full\_name: Klf10

Note: This product is for in vitro research use only