

## Leukotriene B4 receptor 1

Catalog No: #AP84950



Package Size: #AP84950-1 50ug #AP84950-2 100ug #AP84950-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Leukotriene B4 receptor 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LTB4-R 1, Chemoattractant receptor-like 1, G-protein coupled receptor 16, P2Y purinoceptor 7, P2Y7, BLT, BLT1, BLTR, CMKRL1, GPR16, P2RY7
Accession No.	Uniprot ID:Q15722
Uniprot	Q15722
GeneID	1241;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: LTB4R

Note: This product is for in vitro research use only