

Latexin

Catalog No: #AP84951



Package Size: #AP84951-1 50ug #AP84951-2 100ug #AP84951-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Latexin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Endogenous carboxypeptidase inhibitor,ECl,Tissue carboxypeptidase inhibitor,TCI
Accession No.	Uniprot ID:P70202
Uniprot	P70202
GeneID	17035;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lxn

Note: This product is for in vitro research use only