

## Protein LYRIC

Catalog No: #AP84952



Package Size: #AP84952-1 50ug #AP84952-2 100ug #AP84952-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Protein LYRIC
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	3D3/LYRIC,Lysine-rich CEACAM1 co-isolated protein,Metadherin,Metastasis adhesion protein,Lytic
Accession No.	Uniprot ID:Q80WJ7
Uniprot	Q80WJ7
GeneID	67154;
SDS-PAGE MW	27.5
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Expression Region:330-579aaStorage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Mtdh

Note: This product is for in vitro research use only