

Interleukin-17B

Catalog No: #AP84953



Package Size: #AP84953-1 50ug #AP84953-2 100ug #AP84953-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interleukin-17B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	IL-17B,Cytokine CX1,Cytokine-like protein ZCYTO7,Neuronal interleukin-17-related factor,Nirf,Zcyto7
Accession No.	Uniprot ID:Q9QXT6
Uniprot	Q9QXT6
GeneID	56069;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Il17b

Note: This product is for in vitro research use only