

L-lactate dehydrogenase B chain

Catalog No: #AP84956



Package Size: #AP84956-1 50ug #AP84956-2 100ug #AP84956-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	L-lactate dehydrogenase B chain
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LDH-B,LDH heart subunit,LDH-H,Ldh-2,Ldh2
Accession No.	Uniprot ID:P16125
Uniprot	P16125
GeneID	16832;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Ldhb

Note: This product is for in vitro research use only