

Leptin receptor

Catalog No: #AP84957



Package Size: #AP84957-1 50ug #AP84957-2 100ug #AP84957-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leptin receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LEP-R,B219,OB receptor,OB-R,CD295,Db,Obr
Accession No.	Uniprot ID:P48356
Uniprot	P48356
GeneID	16847;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH2PO4, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lepr

Note: This product is for in vitro research use only