

## Arachidonate 12-lipoxygenase, 12S-type

Catalog No: #AP84958



Package Size: #AP84958-1 50ug #AP84958-2 100ug #AP84958-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Arachidonate 12-lipoxygenase, 12S-type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	12S-LOX,Lipoxin synthase 12-LO,Platelet-type lipoxygenase 12,12LO,LOG12
Accession No.	Uniprot ID:P18054
Uniprot	P18054
GeneID	239;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:ALOX12

Note: This product is for in vitro research use only