

Low-density lipoprotein receptor-related protein 8

Catalog No: #AP84962



Package Size: #AP84962-1 50ug #AP84962-2 100ug #AP84962-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor-related protein 8
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Chicken
Other Names	LRP-8, Apolipoprotein E receptor 2, Protein LR8B, LR8B
Accession No.	Uniprot ID:Q98931
Uniprot	Q98931
GeneID	396102;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: LRP8

Note: This product is for in vitro research use only