

Interferon-induced protein 44-like

Catalog No: #AP84965



Package Size: #AP84965-1 50ug #AP84965-2 100ug #AP84965-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Interferon-induced protein 44-like
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	C1orf29,GS3686
Accession No.	Uniprot ID:Q53G44
Uniprot	Q53G44
GeneID	10964;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:IFI44L

Note: This product is for in vitro research use only