

# 6-phosphofructokinase, muscle type

Catalog No: #AP84969

Package Size: #AP84969-1 50ug #AP84969-2 100ug #AP84969-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	6-phosphofructokinase, muscle type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Phosphofructo-1-kinase isozyme A, PFK-A, Phosphofructokinase-M, Phosphofructokinase 1, Phosphohexokinase, Pfk-m, Pfka
Accession No.	Uniprot ID:P47857
Uniprot	P47857
GeneID	18642;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: PfkM

Note: This product is for in vitro research use only