

Mitofusin-2

Catalog No: #AP84970



Package Size: #AP84970-1 50ug #AP84970-2 100ug #AP84970-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mitofusin-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Mitochondrial transmembrane GTPase FZO1A,Protein HSG,Transmembrane GTPase MFN2,Fzo1a
Accession No.	Uniprot ID:Q8R500
Uniprot	Q8R500
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Mfn2

Note: This product is for in vitro research use only