

Keratin, type I cytoskeletal 10

Catalog No: #AP84976



Package Size: #AP84976-1 50ug #AP84976-2 100ug #AP84976-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|--------------------|---|
| Product Name | Keratin, type I cytoskeletal 10 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Dog |
| Other Names | Cytokeratin-10,CK-10,Epithelial keratin-10,Keratin-10,K10,KER10 |
| Accession No. | Uniprot ID:Q6EIZ0 |
| Uniprot | Q6EIZ0 |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT10

Note: This product is for in vitro research use only