

Keratin, type II cytoskeletal 5

Catalog No: #AP84977



Package Size: #AP84977-1 50ug #AP84977-2 100ug #AP84977-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Keratin, type II cytoskeletal 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	58 kDa cytokeratin,Cytokeratin-5,CK-5,Keratin-5,K5,Type-II keratin Kb5
Accession No.	Uniprot ID:P13647
Uniprot	P13647
GeneID	3852;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:KRT5

Note: This product is for in vitro research use only