

Serine protease inhibitor Kazal-type 1

Catalog No: #AP84979



Package Size: #AP84979-1 50ug #AP84979-2 100ug #AP84979-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Serine protease inhibitor Kazal-type 1
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Cholecystokinin-releasing peptide,Monitor peptide,Pancreatic secretory trypsin inhibitor,PSTI-I
Accession No.	Uniprot ID:P09655
Uniprot	P09655
GeneID	24833;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spink1

Note: This product is for in vitro research use only