

Integrin alpha-L

Catalog No: #AP84980



Package Size: #AP84980-1 50ug #AP84980-2 100ug #AP84980-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Integrin alpha-L
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Bovine
Other Names	CD11 antigen-like family member A,Leukocyte adhesion glycoprotein LFA-1 alpha chain,LFA-1A,Leukocyte function-associated molecule 1 alpha chain,CD11a
Accession No.	Uniprot ID:P61625
Uniprot	P61625
GeneID	281874;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:ITGAL

Note: This product is for in vitro research use only