

Transcription factor 4

Catalog No: #AP84981

Package Size: #AP84981-1 50ug #AP84981-2 100ug #AP84981-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|--------------------|--|
| Product Name | Transcription factor 4 |
| Brief Description | Recombinant Protein |
| Host Species | E.coli |
| Purification | Greater than 90% by SDS-PAGE |
| Species Reactivity | Rt |
| Other Names | TCF-4, Immunoglobulin transcription factor 2, ITF-2, R8f DNA-binding protein, SL3-3 enhancer factor 2, SEF-2, Itf2, Sef2 |
| Accession No. | Uniprot ID:Q62655 |
| Uniprot | Q62655 |
| GeneID | 84382; |
| Tag Info | His |
| Storage | Store at -20°C. (Avoid repeated freezing and thawing.) |

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tcf4

Note: This product is for in vitro research use only