

MHC class I polypeptide-related sequence B

Catalog No: #AP84983



Package Size: #AP84983-1 50ug #AP84983-2 100ug #AP84983-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MHC class I polypeptide-related sequence B
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MIC-B,PERB11.2
Accession No.	Uniprot ID:Q29980
Uniprot	Q29980
GeneID	4277;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MICB

Note: This product is for in vitro research use only