

Leukocyte immunoglobulin-like receptor subfamily A member 2

Catalog No: #AP84986

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP84986-1 50ug #AP84986-2 100ug #AP84986-3 1mg

Description

Product Name	Leukocyte immunoglobulin-like receptor subfamily A member 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	CD85 antigen-like family member H,Immunoglobulin-like transcript 1,ILT-1,Leukocyte immunoglobulin-like receptor 7,LIR-7,CD85h,ILT1,LIR7
Accession No.	Uniprot ID:Q8N149
Uniprot	Q8N149
GeneID	11027;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LILRA2

Note: This product is for in vitro research use only