

# Lysophosphatidic acid receptor 3

Catalog No: #AP84987

Package Size: #AP84987-1 50ug #AP84987-2 100ug #AP84987-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Lysophosphatidic acid receptor 3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	LPA receptor 3,Lysophosphatidic acid receptor Edg-7,snGPCR32,Edg7,Lpa3
Accession No.	Uniprot ID:Q8K5E0
Uniprot	Q8K5E0
GeneID	66025;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Lpar3

Note: This product is for in vitro research use only