

Laminin subunit alpha-3

Catalog No: #AP84992



Package Size: #AP84992-1 50ug #AP84992-2 100ug #AP84992-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Laminin subunit alpha-3
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Epiligrin 170 kDa subunit,E170,Epiligrin subunit alpha,Kalinin subunit alpha,Laminin-5 subunit alpha,Laminin-6 subunit alpha,Laminin-7 subunit alpha,Nicein subunit alpha,LAMNA
Accession No.	Uniprot ID:Q16787
Uniprot	Q16787
GeneID	3909;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAMA3

Note: This product is for in vitro research use only