

La-related protein 7

Catalog No: #AP84993

Package Size: #AP84993-1 50ug #AP84993-2 100ug #AP84993-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	La-related protein 7
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	La ribonucleoprotein domain family member 7
Accession No.	Uniprot ID:Q5XI01
Uniprot	Q5XI01
GeneID	686883;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Larp7

Note: This product is for in vitro research use only