

Glycoprotein endo-alpha-1,2-mannosidase

Catalog No: #AP84994



Package Size: #AP84994-1 50ug #AP84994-2 100ug #AP84994-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Glycoprotein endo-alpha-1,2-mannosidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	Endo-alpha mannosidase,Endo-alpha-D-mannosidase,Enman
Accession No.	Uniprot ID:Q5GF25
Uniprot	Q5GF25
GeneID	140808;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Manea

Note: This product is for in vitro research use only