

Ghrelin O-acyltransferase

Catalog No: #AP84995



Package Size: #AP84995-1 50ug #AP84995-2 100ug #AP84995-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ghrelin O-acyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Membrane-bound O-acyltransferase domain-containing protein 4, O-acyltransferase domain-containing protein 4, GOAT, OACT4, FKSG89
Accession No.	Uniprot ID:Q96T53
Uniprot	Q96T53
GeneID	619373;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MBOAT4

Note: This product is for in vitro research use only