

# Intelectin-1a

Catalog No: #AP84996



Package Size: #AP84996-1 50ug #AP84996-2 100ug #AP84996-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Intelectin-1a
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Galactofuranose-binding lectin, Intestinal lactoferrin receptor, Intl, Itln, Itln1a, Itlna
Accession No.	Uniprot ID:O88310
Uniprot	O88310
GeneID	16429;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer: 50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene\_full\_name: Itln1

Note: This product is for in vitro research use only