

Lamina-associated polypeptide 2, isoforms beta/delta/epsilon/gamma

Catalog No: #AP84997

Package Size: #AP84997-1 50ug #AP84997-2 100ug #AP84997-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lamina-associated polypeptide 2, isoforms beta/delta/epsilon/gamma
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Thymopoietin isoforms beta/delta/epsilon/gamma, TP beta/delta/epsilon/gamma, Lap2
Accession No.	Uniprot ID:Q61029
Uniprot	Q61033
GeneID	21917;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Tmpo

Note: This product is for in vitro research use only