

Lymphocyte cytosolic protein 2

Catalog No: #AP84999



Package Size: #AP84999-1 50ug #AP84999-2 100ug #AP84999-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lymphocyte cytosolic protein 2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	SH2 domain-containing leukocyte protein of 76 kDa, SLP-76 tyrosine phosphoprotein, SLP76
Accession No.	Uniprot ID:Q60787
Uniprot	Q60787
GeneID	16822;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Lcp2

Note: This product is for in vitro research use only