

Leukosialin

Catalog No: #AP85001



Package Size: #AP85001-1 50ug #AP85001-2 100ug #AP85001-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Leukosialin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	B-cell differentiation antigen LP-3,Leukocyte sialoglycoprotein,Lymphocyte antigen 48,Ly-48,Sialophorin,CD43
Accession No.	Uniprot ID:P15702
Uniprot	P15702
GeneID	20737;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Spn

Note: This product is for in vitro research use only