

Low-density lipoprotein receptor-related protein 4

Catalog No: #AP85008



Package Size: #AP85008-1 50ug #AP85008-2 100ug #AP85008-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Low-density lipoprotein receptor-related protein 4
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	LRP-4,LDLR dan,Kiaa0816
Accession No.	Uniprot ID:Q8VI56
Uniprot	Q8VI56
GeneID	228357;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lrp4

Note: This product is for in vitro research use only