

HDAC7A (Phospho-Ser486) Antibody

Catalog No: #12685



Package Size: #12685-1 50ul #12685-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HDAC7A (Phospho-Ser486) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WBoO IHC
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of HDAC7A (Phospho Ser486) Rabbit pAb at Human, Mouse,Rat
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HDAC7A (Phospho-Ser486)
Target Name	HDAC7A
Modification	Phospho
Other Names	HDAC7A (Phospho-Ser486) Antibody
Accession No.	Swiss-Prot#: Q8WUI4NCBI Gene ID: 51564
Uniprot	Q8WUI4
GeneID	51564;
Target Species	human
Calculated MW	103kD
Concentration	1.0mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	-15°C to -25°C/1 year(Do not lower than -25°C)

Application Details

WB 1:500-2000 IHC 1:50-200

Note: This product is for in vitro research use only