

Protein kinase C epsilon type

Catalog No: #AP85010

Package Size: #AP85010-1 50ug #AP85010-2 100ug #AP85010-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Protein kinase C epsilon type
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	nPKC-epsilon,PKCE
Accession No.	Uniprot ID:Q02156
Uniprot	Q02156
GeneID	5581;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PRKCE

Note: This product is for in vitro research use only