

Lymphocyte activation gene 3 protein

Catalog No: #AP85011



Package Size: #AP85011-1 50ug #AP85011-2 100ug #AP85011-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lymphocyte activation gene 3 protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	LAG-3,Protein FDC,CD223,FDC
Accession No.	Uniprot ID:P18627
Uniprot	P18627
GeneID	3902;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:LAG3

Note: This product is for in vitro research use only