

## Laminin subunit beta-2

Catalog No: #AP85012



Package Size: #AP85012-1 50ug #AP85012-2 100ug #AP85012-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Laminin subunit beta-2
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Laminin B1s chain,Laminin-11 subunit beta,Laminin-14 subunit beta,Laminin-15 subunit beta,Laminin-3 subunit beta,Laminin-4 subunit beta,Laminin-7 subunit beta,Laminin-9 subunit beta,S-laminin subunit beta,S-LAM beta,LAMS
Accession No.	Uniprot ID:P55268
Uniprot	P55268
GeneID	3913;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:LAMB2

Note: This product is for in vitro research use only