

Lutropin subunit beta

Catalog No: #AP85014



Package Size: #AP85014-1 50ug #AP85014-2 100ug #AP85014-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Lutropin subunit beta
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lutropin beta chain,Luteinizing hormone subunit beta,LH-B
Accession No.	Uniprot ID:O09108
Uniprot	O09108
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lhb

Note: This product is for in vitro research use only