

Ghrelin O-acyltransferase

Catalog No: #AP85019



Package Size: #AP85019-1 50ug #AP85019-2 100ug #AP85019-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ghrelin O-acyltransferase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Membrane-bound O-acyltransferase domain-containing protein 4, O-acyltransferase domain-containing protein 4, Gm171, Goat
Accession No.	Uniprot ID:P0C7A3
Uniprot	P0C7A3
GeneID	234155;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Mboat4

Note: This product is for in vitro research use only