

## Mitochondrial carnitine/acylcarnitine carrier protein

Catalog No: #AP85020



Package Size: #AP85020-1 50ug #AP85020-2 100ug #AP85020-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mitochondrial carnitine/acylcarnitine carrier protein
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Carnitine/acylcarnitine translocase,CAC,Solute carrier family 25 member 20,CAC,CACT
Accession No.	Uniprot ID:O43772
Uniprot	O43772
GeneID	788;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SLC25A20

Note: This product is for in vitro research use only