

## Hepatocyte growth factor receptor

Catalog No: #AP85021



Package Size: #AP85021-1 50ug #AP85021-2 100ug #AP85021-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Hepatocyte growth factor receptor
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Rt
Other Names	HGF receptor,HGF/SF receptor,Proto-oncogene c-Met,Scatter factor receptor,SF receptor,Tyrosine-protein kinase Met
Accession No.	Uniprot ID:P97523
Uniprot	P97523
GeneID	24553;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:Met

Note: This product is for in vitro research use only