

Microfibrillar-associated protein 5

Catalog No: #AP85022



Package Size: #AP85022-1 50ug #AP85022-2 100ug #AP85022-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Microfibrillar-associated protein 5
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	MFAP-5,MP25, Microfibril-associated glycoprotein 2, MAGP-2, MAGP2
Accession No.	Uniprot ID:Q13361
Uniprot	Q13361
GeneID	8076;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: MFAP5

Note: This product is for in vitro research use only