

Pancreatic triacylglycerol lipase

Catalog No: #AP85025



Package Size: #AP85025-1 50ug #AP85025-2 100ug #AP85025-3 1mg

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Pancreatic triacylglycerol lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	PL
Accession No.	Uniprot ID:P16233
Uniprot	P16233
GeneID	5406;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:PNLIP

Note: This product is for in vitro research use only