

Pancreatic triacylglycerol lipase

Catalog No: #AP85026



Package Size: #AP85026-1 50ug #AP85026-2 100ug #AP85026-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Pancreatic triacylglycerol lipase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	PL
Accession No.	Uniprot ID:Q6P8U6
Uniprot	Q6P8U6
GeneID	69060;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer: 50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole, 10% glycerol (pH 8.0) gene_full_name: Pnlip

Note: This product is for in vitro research use only