

Protein-lysine 6-oxidase

Catalog No: #AP85027



Package Size: #AP85027-1 50ug #AP85027-2 100ug #AP85027-3 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Protein-lysine 6-oxidase
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Ms
Other Names	Lysyl oxidase,Ras excision protein,Rrg
Accession No.	Uniprot ID:P28301
Uniprot	P28301
GeneID	16948;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

Application Details

Storage Buffer:50mM NaH₂PO₄, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene_full_name:Lox

Note: This product is for in vitro research use only