

# Kallistatin

Catalog No: #AP85028



Package Size: #AP85028-1 50ug #AP85028-2 100ug #AP85028-3 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Kallistatin
Brief Description	Recombinant Protein
Host Species	E.coli
Purification	Greater than 90% by SDS-PAGE
Species Reactivity	Hu
Other Names	Kallikrein inhibitor,Peptidase inhibitor 4,PI-4,Serpin A4,KST,PI4
Accession No.	Uniprot ID:P29622
Uniprot	P29622
GeneID	5267;
Tag Info	His
Storage	Store at -20°C. (Avoid repeated freezing and thawing.)

## Application Details

Storage Buffer:50mM NaH<sub>2</sub>PO<sub>4</sub>, 500mM NaCl Buffer with 500mM Imidazole,10%glycerol(PH8.0)gene\_full\_name:SERPINA4

Note: This product is for in vitro research use only